



#### **CERVICAL CANCER NOW HAS SOLUTION FOR PREVENTION**

#### Screen & Treat in a single sitting!

Elimination of cervical cancer by 2030



In 2019, WHO issued new guidelines for the use of thermal ablation for cervical pre-cancerous lesions, recommending the thermal ablation as an alternative method of treatment.

Thermo Glide

Treating precancerous lesions

with thermocoagulation

**TREAT** 



#### **SCREEN**
























### **EVA PRO-**

#### The Gynac Stethoscope

A turn-keycervicaland vulvar visualization solution to streamline how you perform, documentand review patient cases.

#### Digital touch screen

Enhances the visualization experience, allowing to zoom in on images and perform the examination in a convenient manner.



Save images with notes, annotate biopsies and points of interest, add clinical impressions, edit in device or on portal.

#### Built in camera with magnifying

you can identify a lesion as small as 0.05 mm



Hands free image capture For optimal use









## capabilities

Superior image quality,

#### Powered by Al



supporting the WHO strategic objective to eliminate cervical cancer by 2030



can be used by trained doctors & paramedics



Results received within 60 seconds, directly at point of care





## Pre-cancerous lesions identified by EVA PRO can now be treated on site

# ThermoGlide

# An immediate point-of-care treatment option to prevent cervical cancer

## A complete thermo coagulation solution for the treatment of precancerous cervical lesions

The ThermoGlide uses heat to decimate tissue. The probe is heated quickly within **8 seconds and remains at 100°C** allowing for rapid cervical ablation. The superficial epithelium withers and dissipates after the treatment, while the stroma and glandular crypts are destroyed.

Thermo coagulation offers shorter treatment time as compared to cryotherapy.



- 4 probes
  Versatile treatment entire
  - Versatile treatment options
- ② 2 battery packs
  Wireless battery pack
- Charger

  Rechargeable battery, with no need for electricity
- Minimal pain and discomfort
  Compared to the conventional treatment methods
- Safe and effective
   Timer function & constant temperature
   setting reducing possibility of human error
- Quick learning curve, simple handling and assembly



In 2019, WHO issued new guidelines for the use of thermal ablation for cervical pre-cancerous lesions, recommending the thermal ablation as an alternative method of treatment.





# SPECIAL OFFERS NOW AVAILABLE

# BUYBACK EXISTING COLPOSCOPE TO EVAPro

Al Enabled Digital Colposcope

# INTRODUCING ThermoGlide

An immediate point-of-care treatment option precancerous lesions with thermocoagulation





- \* CIDEX VOUCHER WORTH ₹8,500 FREE
- \* AVAIL ATTRACTIVE BUY BACK OFFERS
- **\* 0 DOWN PAYMENT**

Pay in equal instalments of 12 months, 18 months, or 24 months!

SCAN TO KNOW MORE





