

MAXIMIZE YOUR CHANCE OF FIRST-ATTEMPT SUCCESS.¹

You prepare for every eventuality before you intubate.

So when your patient is on the table, when those critical moments are passing, you are better prepared for whatever happens.

The McGRATH[®] MAC enhanced direct laryngoscope (EDL) is a single handheld device that combines the benefits of:

- Direct laryngoscopy (DL) — the familiar technique of a Macintosh laryngoscope
- Video laryngoscopy (VL) — video technology

THE McGRATH[®] MAC ENHANCED DIRECT LARYNGOSCOPE

INTUITIVE

Uses standard DL technique combined with VL technology

PATIENT-FRIENDLY

Offers improved access to the patient's airway with a slim blade profile²

INSTILLS CONFIDENCE

Includes a 250-minute battery with a minute-by-minute on-screen charge countdown that confirms device battery life

DURABLE

Drop tested to two meters, the McGRATH[®] MAC EDL can endure rough handling and still perform in a fast-paced medical environment

EASY TO CLEAN

Provides sterile integrity with disposable blades and is fully immersible. Both screen and handle can be treated with high-level disinfection



ORDERING INFORMATION

SKU NO.	DESCRIPTION	SIZE	QUANTITY
300-000-000	Handle	N/A	1
350-017-000	McGRATH [®] MAC blade size 2	2	50/box
350-005-000	McGRATH [®] MAC blade size 3	3	50/box
350-013-000	McGRATH [®] MAC blade size 4	4	50/box
X3-003-000	McGRATH [®] MAC X blade size 3	3	10/box
340-000-000	250-minute battery pack	N/A	1

TECHNICAL SPECIFICATIONS

McGRATH[®] MAC Handle

Size: 180 mm x 68 mm x 110 mm

Weight: 200 g

Power: Proprietary 3.6V lithium battery pack (c. 250 minutes)

Light source: High-intensity LED

Display: 2.5" LCD color display

Camera: CMOS


Materials: Durable medical-grade thermoplastics with reinforced structural alloy core. The device and packaging are latex-free.

McGRATH[®] MAC Disposable Blades

Material: Fog-free medical-grade optical polymer

Packaging: Packaged sterile for single use

For demonstrations and quotes, please contact your local representative. Find out more at www.medtronic.com/covidien/products





INTUBATION IS EASY. UNTIL IT ISN'T.

McGRATH[®] MAC
Enhanced Direct Laryngoscope

Medtronic
Further, Together



1. Su YC, Chen CC, Lee YK, Lee JY, Lin KJ. Comparison of video laryngoscopes with direct laryngoscopy for tracheal intubation: a meta-analysis of randomised trials. *Eur J Anaesthesiol.* 2011;28(11):788-795.
2. Levitan R, Heitz J, Sweeney M, Cooper R. The complexities of tracheal intubation with direct laryngoscopy and alternative intubation devices. *Ann Emerg Med.* 2011;57 (3):240-247.

©2013, 2015, Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. TM® Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

Medtronic

